

DocuMentor: Interactive PDF Knowledge Bot

Upload your PDF file



Drag and drop file here

Limit 200MB per file • PDF

Browse files



Please upload a PDF to start interacting with your smart chatbot.

DocuMentor: Interactive PDF Knowledge Bot

Upload your PDF file



Drag and drop file here

Limit 200MB per file • PDF

Browse files



H-046-029216-00-4.0-hybase-v6-user-manual-en.pdf 4.7MB



Hang Tight Smarty! Processing...

>>

Stop Deploy

What is the procedure for disinfecting or cleaning the table?

The cleaning procedure for the operating table is detailed on page 93. Before cleaning, you must adjust the table to the 0 position, lift the tabletop to its highest position, and switch off and disconnect the power cord [page 93].

The cleaning process itself involves these steps [page 93]:

1. Gently and thoroughly wipe off dust using a lint-free cloth.
2. Wipe the surface with a lint-free cloth moistened with a cleaner (clean water or a soap solution).
3. Clean the surface again with a disposable, lint-free cloth moistened with water.
4. Finally, wipe the surface dry with a dry lint-free cloth.

Important considerations from page 92 remind us to:

- **Avoid alcohol-based cleaners and disinfectants** near high-frequency surgical equipment.
- **Always wear gloves, masks, and eye protection.**
- **Never put your hands into gaps** in the tabletop during cleaning or disinfection.
- **Refer to the instructions for use** of your specific cleaning agents and disinfectants. Cleaning and disinfection of other reusable accessories should also follow their specific instructions.

Remember that the responsible hospital or institution is ultimately responsible for all cleaning and disinfection procedures [page 92].

Ask me something about the uploaded PDF...



>>

Deploy

Important considerations from page 92 remind us to:

- **Avoid alcohol-based cleaners and disinfectants** near high-frequency surgical equipment.
- **Always wear gloves, masks, and eye protection.**
- **Never put your hands into gaps** in the tabletop during cleaning or disinfection.
- **Refer to the instructions for use** of your specific cleaning agents and disinfectants. Cleaning and disinfection of other reusable accessories should also follow their specific instructions.

Remember that the responsible hospital or institution is ultimately responsible for all cleaning and disinfection procedures [page 92].

> Sources

Why to avoid alcohol?

Avoid alcohol-based cleaners and disinfectants because they can create explosive vapor mixtures and ignite near high-frequency surgical equipment [page 92, page 13]. This poses a significant fire hazard.

▼ Sources

- Page 13 → cause personal injury. If accessories of other manufacturers must be used, their use should be expressly permitted by Mindray. • Before using accessories, be sure to read the operator's manual thorou...

Ask me something about the uploaded PDF...



Deploy

Session ID: c27bffb6-0632-4b36-b5b6-885535992d11

Full Chat History

History

User: What is the procedure for disinfecting or cleaning the table?

Bot: The cleaning procedure for the operating table is detailed on page 93. Before cleaning, you must adjust the table to the 0 position, lift the tabletop to its highest position, and switch off and disconnect the power cord [page 93].

The cleaning process itself involves these steps [page 93]:

1. Gently and thoroughly wipe off dust using a lint-free cloth.
2. Wipe the surface with a lint-free cloth moistened with a cleaner (clean water or a soap solution).
3. Clean the surface again with a disposable, lint-free cloth moistened with water.
4. Finally, wipe the surface dry with a dry lint-free cloth.

Important considerations from page 92 remind us to:

- Page 13 → cause personal injury. If accessories of other manufacturers must be used, their use should be expressly permitted by Mindray. • Before using accessories, be sure to read the operator's manual thorou...
- Page 92 → 5 - 1 5 Cleaning and Disinfection 5.1 Safety Information The product must be cleaned and disinfected after every use. WARNING • Use only Mindray approved cleaners and disinfectants and methods listed ...



tell me the page number for the sterilization procedure.



The provided text mentions sterilization on page 108, but it only states that the operating table is *not* sterilizable and gives operational limitations. There is no sterilization *procedure* detailed on that page or in the provided text.

Sources

- Page 108 → of sterilization Non-sterilization. According to the suitability for use in an oxygen rich environment The operating table cannot be used in an oxygen rich environment. According to the mode of op...

Ask me something about the uploaded PDF...

